

EUROPEAN AND MEDITERRANEAN PLANT PROTECTION ORGANIZATION
ORGANISATION EUROPEENNE ET MEDITERRANEENNE POUR LA PROTECTION DES PLANTES
Summary sheet of validation data for a diagnostic test

The EPPO Standard PM 7/98 *Specific requirements for laboratories preparing accreditation for a plant pest diagnostic activity* describes how validation should be conducted. It also includes definitions of performance criteria.

Laboratory contact details	Anses Plant Health Laboratory - Pests and Tropical Pathogens Unit Pôle de Protection des Plantes, 7 Chemin de l'IRAT, 97410 Saint Pierre, France
Short description of the test	Detection of <i>Xanthomonas axonopodis</i> pv. <i>dieffenbachiae</i> by DAS-ELISA in leaves and pure culture
Date, reference of the validation report	2012-03-01 - Inter-laboratory ring test : <i>Xanthomonas axonopodis</i> pv. <i>dieffenbachiae</i> in Anthurium (Report Xad01-version 2)
Validation process according to EPPO Standard PM7/98?	no
Is the lab accredited for this test?	no
Was the validated data generated in the framework of a project?	
Description of the test	
Organism(s)	<i>Xanthomonas axonopodis</i> pv. <i>dieffenbachiae</i> (XANTDF)
Detection / identification	detection
Method(s)	Extraction Serological DAS-ELISA
Method: Extraction	
Reference of the test description	
As or adapted from an EPPO diagnostic protocol	yes
EPPO Diagnostic Protocol name	PM 7/023 <i>Xanthomonas axonopodis</i> pv. <i>dieffenbachiae</i> (version 2)
Name of the test	Extraction from symptomatic plant material in PBS buffer (Appendix 1.1)
Other information	
Other details on the test	Extraction as in Appendix 1 of PM7/23(2)
Method: Serological DAS-ELISA	
Reference of the test description	
As or adapted from an EPPO diagnostic	yes

protocol	
EPPO Diagnostic Protocol name	PM 7/023 Xanthomonas axonopodis pv. dieffenbachiae (version 2)
Name of the test	DAS-ELISA
<i>Other information</i>	
Other details on the test	PRI protocol for detection of Xad by DAS-ELISA
Are the performance characteristics included in the EPPO diagnostic protocol?	no
Performance Criteria :	
Organism 1.:	Xanthomonas axonopodis pv. dieffenbachiae(XANTDF)
<u>Analytical sensitivity</u>	
What is smallest amount of target that can be detected reliably?	10 ⁴ CFU.mL-1
<u>Diagnostic sensitivity</u>	
Proportion of infected/infested samples tested positive compared to results from the standard test, see appendix 2 of PM 7/98	Comparative study : 100% ; Collaborative study : 63.5%-70.5%
Standard test(s)	Isolation + AGDIA Indirect-ELISA on pure culture (OEPP PM7/23)
<u>Analytical specificity - inclusivity</u>	
Number of strains/populations of target organisms tested	50 (see attached downloadable file Appendix 1)
Specificity value	100
<u>Analytical specificity - exclusivity</u>	
Number of non-target organisms tested	53 (see attached downloadable file Appendix 2)
Specificity value	exclusivity = 0,47. DAS-ELISA can not exclude strains described as Xad but not pathogenic to Anthurium and some strains that belong to the same species but to a different pathovar. In addition, this method can not exclude some saprophytic strains.
<u>Diagnostic Specificity</u>	
Proportion of uninfected/uninfested samples (true negatives) testing negative compared to results from a standard test	Comparative study : 94%; Collaborative study : 91%-93%
Specify the test(s)	Isolation + AGDIA Indirect-ELISA on pure culture (OEPP PM7/23)
<u>Reproducibility</u>	
Provide the calculated % of agreement for a given level of the pest (see PM 7/98)	74%-78%
<u>Repeatability</u>	
Provide the calculated % of agreement for a given level of the pest (see PM 7/98)	89%-91%

Test performance study	
Test performance study?	yes
Brief details of the test performance study and its output. It is available, link to published article/report	Because of its lack of specificity, the DAS-ELISA should not be used for identification purposes. Report is downloadable from this database.
The following complementary files are available online:	<ul style="list-style-type: none"> • Appendix 1-List target strains-2 • Appendix 2-List non target strains-2 • EILVReport-V02 01.03.2012 correction

Creation date: 2012-05-14 00:00:00 - Last update: 2020-09-14 12:31:26