

**EUROPEAN AND MEDITERRANEAN PLANT PROTECTION ORGANIZATION**  
**ORGANISATION EUROPEENNE ET MEDITERRANEENNE POUR LA PROTECTION DES PLANTES**  
**Summary sheet of validation data for a diagnostic test**

The EPPO Standard PM 7/98 *Specific requirements for laboratories preparing accreditation for a plant pest diagnostic activity* describes how validation should be conducted. It also includes definitions of performance criteria.

<b>Laboratory contact details</b>	Bacteriology. Instituto Valenciano de Investigaciones Agrarias CV-315, km. 10.7, 46113 Moncada, Spain
<b>Short description of the test</b>	Proficiency testing for the evaluation of molecular diagnosis of <i>Xylella fastidiosa</i> . The panel of samples included randomized <i>Xylella</i> -free preparations and samples spiked with heat-inactivated (incubation at 70C for 15min) bacterial suspensions at different concentrations, prepared by scraping 10 days old colonies of <i>X. fastidiosa</i> subsp. <i>pauca</i> strain CoDiRO (CFBP 8402)
<b>Date, reference of the validation report</b>	2021-03-08 - EU-XF- PT-2017-02
<b>Validation process according to EPPO Standard PM7/98?</b>	yes
<b>Is the lab accredited for this test?</b>	no
<b>Was the validated data generated in the framework of a project?</b>	Other_project
<b>If yes, please specify</b>	POnte project
<b>Description of the test</b>	
<b>Organism(s)</b>	<i>Xylella fastidiosa</i> (XYLEFA)
<b>Detection / identification</b>	detection
<b>Method(s)</b>	Molecular Extraction DNA RNA Molecular Conventional PCR Molecular real time PCR
<b>Method: Molecular Extraction DNA RNA</b>	
<i>Reference of the test description</i>	
<i>Other information</i>	
<b>Method: Molecular Conventional PCR</b>	
<i>Reference of the test description</i>	
<i>Other information</i>	
<b>Method: Molecular real time PCR</b>	
<i>Reference of the test description</i>	
<i>Other information</i>	
<b>Performance Criteria :</b>	

<b>Organism 1.:</b>	<b>Xylella fastidiosa(XYLEFA)</b>
<b>Test performance study</b>	
<b>Test performance study?</b>	no
The following complementary files are available online:	<ul style="list-style-type: none"> <li>• <a href="#">Proficiency testing for the evaluation of molecular and serological diagnosis of Xylella fastidiosa</a></li> </ul>

Creation date: 2021-03-08 09:10:02 - Last update: 2021-03-08 09:13:19